

# Simulation of the Lungs for the Respiration Process

Bharath Mathaji Chellaboyina

International Institute of Information Technology


Rajiv Gandhi University of Knowledge Technologies, Andhra Pradesh

## Background:

Perhaps the most difficult aspect of COVID-19 is that a recently contaminated individual shows diagnosable manifestations, for example, internal heat level ascent, a few days subsequent to getting the infection. In the beginning stage of disease (i.e., hatching or incubation period), an undiscovered and uninformed individual can spread the infection to other people. The quickest and most effective course of Coronavirus transmission is the respiratory framework. In this manner, fostering a model of the respiratory framework to foresee changes in the lung execution upon Coronavirus contamination is helpful for early conclusion and mediation during the hatching time frame.

## Description of the flowsheet:

A simple model was developed by performing mass and energy balances on the lungs; it was simulated by the DWSIM chemical process simulator using Peng Robinson Method. S-01 (Inhaled Humid Air) Consists of 20<sup>o</sup> C at 101.3 KPa having 0.01607 Kmol/hr molar flow with mole fractions Oxygen of 0.2076, Nitrogen of 0.7808 and water of 0.0116 And E-01(Sensible Heat Exchange) were fed into the HEAT-01 (Heater). The heater has S-03 (Warm Humid Air) with  $\Delta P$  0 Kpa and  $\Delta T$  17<sup>o</sup> C in the heater resulting in the temperature rise to 37<sup>o</sup>C and the pressure remains constant. S-02 (Co<sub>2</sub> from the blood) having 37<sup>o</sup>C at 101.3 Kpa and 0.00142 Kmol/hr with Co<sub>2</sub> mole fraction 1, S-04 (Moisture from the lung tissue) having 37<sup>o</sup>C at 101.3 Kpa and 0.0008862 Kmol/hr with H<sub>2</sub>O mole fraction 1, E-02 (Evaporative Heat Exchange) were fed into SEP-01 (Humidifier) with operation conditions 37<sup>o</sup>C at 101.3 Kpa having S-05 (Gas from the lungs) and S-06 (Liquid from the lungs) as outlet streams. Finally vapour and liquid have been separated. S-05 is fed into a CS-01 (Gas Separator) with S-07 (Exhaled Air) and S-08 (O<sub>2</sub> to blood) used oxygen fraction of 1 in outlet stream 2 as spec.



Simulation of the Lungs for the Respiration Process

**Results:**

Object	S-01	S-02	S-03	S-04	
Temperature	20	37	37	37	C
Pressure	1.013	1.013	1.013	1.013	bar
Molar Flow	0.01607	0.00142	0.01607	0.0008862	kmol/h
Molar Flow (Vapor)	0.01607	0.00142	0.01607	0	kmol/h
Molar Fraction (Vapor) / Nitrogen	0.7808	0	0.7808	0	
Molar Fraction (Vapor) / Oxygen	0.2076	0	0.2076	0	
Molar Fraction (Vapor) / Carbon dioxide	0	1	0	0	
Molar Fraction (Vapor) / Water	0.0116	0	0.0116	1	
Object	S-05	S-06	S-07	S-08	
Temperature	37	37	37	37	C
Pressure	1.013	1.013	1.013	1.013	bar
Molar Flow	0.0182977	7.85E-05	0.0182644	3.34E-05	kmol/h
Molar Flow (Vapor)	0.0182977	0	0.0182624	3.34E-05	kmol/h
Molar Fraction (Vapor) / Nitrogen	0.685738	0	0.687065	0	
Molar Fraction (Vapor) / Oxygen	0.182325	0	0.180851	1	
Molar Fraction (Vapor) / Carbon dioxide	0.0776051	0	0.0777552	0	
Molar Fraction (Vapor) / Water	0.0543319	0	0.0543287	0	

**Conclusion:**

To make up for the infection tainted lung failure, the O2 focus expanded in the breathed out air at the expense of diminished CO2 fixation. In opposition to past discoveries on the diminished soundness of Covid in hot and sticky conditions, it was observed that exceptionally hot and moist conditions advance the viral transmission rate in view of the immediate hotness move to the body through breath and buildup of water fumes that might cause disease in the respiratory tract.